



PH-W420

Thermostatic Water Bath Cauldron



PH-W600

Thermostatic Water Bath Cauldron

Technical Parameter

Model	420	600
Heating Power	1000W	1500W
Volume	420*180*200	600*300*220
Power Supply	220V 50HZ	220V 50HZ
Constant Temperature Range	RT~100 °C	RT~100 °C

Features

- This product adopts the sink type, the inner tank is welded with imported stainless steel plates, the outer shell is made of high-quality cold board and sprayed with plastic, and the surface is bright and wear-resistant.
- The electric heating tube is sandwiched in the middle of the water, which has fast heating, high thermal efficiency and low power consumption.
- Digital display temperature control RT ~ 100 °C , according to the needs of use, you can adjust the intelligent temperature controller to control the temperature.